


SAGE[®] PEEK

CERVICAL INTERBODY
FUSION SYSTEM


CERVICAL INTERBODY

Features & Benefits



SAGE® PEEK CERVICAL INTERBODY FUSION SYSTEM

The Sage® Cervical Interbody Fusion System is designed to support cervical loads while maximizing surface contact between the implant and the vertebral bodies. The interbody is offered in four footprints to accommodate individual patient anatomy while the streamlined instrumentation is designed to aid in consistent and accurate implant placement.



System Features:

- Manufactured of biocompatible PEEK polymer
- 12x11, 14x12, 16x13 and 18x14mm footprints in 5-12mm heights, offered in 1mm increments with 7° lordosis
- Tantalum radiopaque markers to optimize visibility and placement
- Multi-directional, aggressive teeth designed to prevent implant migration
- Robust inserter engagement features to facilitate implant insertion
- Large, central graft cavity for packing bone graft to facilitate fusion
- Lateral windows provide vascularization to the fusion site
- Color-coded instruments for easy identification

PEEK Cervical Interbodies

Catalog Number	Description (W x D x H)	Graft Volume	Catalog Number	Description (W x D x H)	Graft Volume
C204-121105-7	12 x 11 x 5mm	.16cc	C204-161305-7	16 x 13 x 5mm	.34cc
C204-121106-7	12 x 11 x 6mm	.20cc	C204-161306-7	16 x 13 x 6mm	.42cc
C204-121107-7	12 x 11 x 7mm	.24cc	C204-161307-7	16 x 13 x 7mm	.50cc
C204-121108-7	12 x 11 x 8mm	.28cc	C204-161308-7	16 x 13 x 8mm	.58cc
C204-121109-7	12 x 11 x 9mm	.32cc	C204-161309-7	16 x 13 x 9mm	.66cc
C204-121110-7	12 x 11 x 10mm	.36cc	C204-161310-7	16 x 13 x 10mm	.74cc
C204-121111-7	12 x 11 x 11mm	.39cc	C204-161311-7	16 x 13 x 11mm	.82cc
C204-121112-7	12 x 11 x 12mm	.43cc	C204-161312-7	16 x 13 x 12mm	.91cc
C204-141205-7	14 x 12 x 5mm	.25cc	C204-181405-7	18 x 14 x 5mm	.44cc
C204-141206-7	14 x 12 x 6mm	.31cc	C204-181406-7	18 x 14 x 6mm	.55cc
C204-141207-7	14 x 12 x 7mm	.37cc	C204-181407-7	18 x 14 x 7mm	.66cc
C204-141208-7	14 x 12 x 8mm	.43cc	C204-181408-7	18 x 14 x 8mm	.77cc
C204-141209-7	14 x 12 x 9mm	.49cc	C204-181409-7	18 x 14 x 9mm	.88cc
C204-141210-7	14 x 12 x 10mm	.55cc	C204-181410-7	18 x 14 x 10mm	.99cc
C204-141211-7	14 x 12 x 11mm	.61cc	C204-181411-7	18 x 14 x 11mm	1.10cc
C204-141212-7	14 x 12 x 12mm	.67cc	C204-181412-7	18 x 14 x 12mm	1.21cc



For more information or to place an order call: 877.9CURITEVA, email: customersupport@curiteva.com or visit www.curiteva.com

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2021 by Curiteva, Inc.
QMSF-1126 Rev 3.0