


SAVANT™ PEEK

POSTERIOR LUMBAR
INTERBODY FUSION SYSTEM


LUMBAR INTERBODY

Features & Benefits



CURITEVA®

SAVANT™ PEEK PLIF SYSTEM



The Savant PEEK Posterior Lumbar Interbody Fusion System is offered in a variety of lengths, heights and lordotic angles providing an ideal fit for every patient. The superior and inferior surfaces of the construct have a pattern of teeth designed to provide increased stability and to help prevent movement of the device. The Savant Interbody Fusion System features instruments that are simple and intuitive to use.

System Features:

- Manufactured out of biocompatible PEEK polymer
- Contoured, distracting nose to facilitate placement
- Central and lateral graft windows for biologic seeding
- Convex design offers a more accurate fit to the patient's anatomy
- Multi-directional teeth resist migration
- Strategically positioned tantalum markers for greater visualization and placement
- Comprehensive disectomy instruments
- Full range of Trials for proper sizing of implant height

PEEK PLIF Interbodies

22mm Length, 10mm Width, 0°		26mm Length, 10mm Width, 0°		30mm Length, 10mm Width, 0°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
6mm	C304-221006-0	6mm	C304-261006-0	6mm	C304-301006-0
7mm	C304-221007-0	7mm	C304-261007-0	7mm	C304-301007-0
8mm	C304-221008-0	8mm	C304-261008-0	8mm	C304-301008-0
9mm	C304-221009-0	9mm	C304-261009-0	9mm	C304-301009-0
10mm	C304-221010-0	10mm	C304-261010-0	10mm	C304-301010-0
11mm	C304-221011-0	11mm	C304-261011-0	11mm	C304-301011-0
12mm	C304-221012-0	12mm	C304-261012-0	12mm	C304-301012-0
13mm	C304-221013-0	13mm	C304-261013-0	13mm	C304-301013-0
14mm	C304-221014-0	14mm	C304-261014-0	14mm	C304-301014-0
15mm	C304-221015-0	15mm	C304-261015-0	15mm	C304-301015-0
16mm	C304-221016-0	16mm	C304-261016-0	16mm	C304-301016-0

22mm Length, 10mm Width, 8°		26mm Length, 10mm Width, 8°		30mm Length, 10mm Width, 8°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
7mm	C304-221007-8	8mm	C304-261008-8	9mm	C304-301009-8
8mm	C304-221008-8	9mm	C304-261009-8	10mm	C304-301010-8
9mm	C304-221009-8	10mm	C304-261010-8	11mm	C304-301011-8
10mm	C304-221010-8	11mm	C304-261011-8	12mm	C304-301012-8
11mm	C304-221011-8	12mm	C304-261012-8	13mm	C304-301013-8
12mm	C304-221012-8	13mm	C304-261013-8	14mm	C304-301014-8
13mm	C304-221013-8	14mm	C304-261014-8	15mm	C304-301015-8
14mm	C304-221014-8	15mm	C304-261015-8	16mm	C304-301016-8
15mm	C304-221015-8	16mm	C304-261016-8		
16mm	C304-221016-8				



For more information or to place an order call: 877.9CURITEVA, email: customersupport@curiteva.com or visit www.curiteva.com

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2020 by Curiteva, Inc.
QMSF-1146 Rev 1.0