


SAVANT™ PEEK

TRANSFORAMINAL LUMBAR
INTERBODY FUSION SYSTEM


LUMBAR INTERBODY

Features & Benefits



CURITEVA®

SAVANT™ PEEK TLIF SYSTEM



The Savant Transforaminal Lumbar Interbody Fusion System anatomical design delivers intraoperative versatility and performance. The interbody is offered in multiple footprints, heights and lordosis providing an ideal fit for every patient. The Savant Interbody Fusion System features instrumentation that allows for better control and more precise implant placement.

System Features:

- Convex design offers a more accurate fit to the patient's anatomy
- Central and lateral graft windows for biologic seeding
- Contoured, distracting nose to facilitate placement
- Multi-directional teeth resist migration
- Angled lateral windows on anterior side of implant to retain bone graft material
- Strategically positioned tantalum markers for greater visualization and placement
- Comprehensive discectomy instruments
- Full range of Trials for proper sizing of implant height

PEEK TLIF Interbodies

28mm Length, 9mm Width, 0°		28mm Length, 11mm Width, 0°		32mm Length, 9mm Width, 0°		32mm Length, 11mm Width, 0°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
6mm	C404-280906-0	6mm	C404-281106-0	6mm	C404-320906-0	6mm	C404-321106-0
7mm	C404-280907-0	7mm	C404-281107-0	7mm	C404-320907-0	7mm	C404-321107-0
8mm	C404-280908-0	8mm	C404-281108-0	8mm	C404-320908-0	8mm	C404-321108-0
9mm	C404-280909-0	9mm	C404-281109-0	9mm	C404-320909-0	9mm	C404-321109-0
10mm	C404-280910-0	10mm	C404-281110-0	10mm	C404-320910-0	10mm	C404-321110-0
11mm	C404-280911-0	11mm	C404-281111-0	11mm	C404-320911-0	11mm	C404-321111-0
12mm	C404-280912-0	12mm	C404-281112-0	12mm	C404-320912-0	12mm	C404-321112-0
13mm	C404-280913-0	13mm	C404-281113-0	13mm	C404-320913-0	13mm	C404-321113-0
14mm	C404-280914-0	14mm	C404-281114-0	14mm	C404-320914-0	14mm	C404-321114-0
15mm	C404-280915-0	15mm	C404-281115-0	15mm	C404-320915-0	15mm	C404-321115-0
16mm	C404-280916-0	16mm	C404-281116-0	16mm	C404-320916-0	16mm	C404-321116-0

28mm Length, 9mm Width, 8°		28mm Length, 11mm Width, 8°		32mm Length, 9mm Width, 8°		32mm Length, 11mm Width, 8°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
7mm	C404-280907-8	7mm	C404-281107-8	7mm	C404-320907-8	7mm	C404-321107-8
8mm	C404-280908-8	8mm	C404-281108-8	8mm	C404-320908-8	8mm	C404-321108-8
9mm	C404-280909-8	9mm	C404-281109-8	9mm	C404-320909-8	9mm	C404-321109-8
10mm	C404-280910-8	10mm	C404-281110-8	10mm	C404-320910-8	10mm	C404-321110-8
11mm	C404-280911-8	11mm	C404-281111-8	11mm	C404-320911-8	11mm	C404-321111-8
12mm	C404-280912-8	12mm	C404-281112-8	12mm	C404-320912-8	12mm	C404-321112-8
13mm	C404-280913-8	13mm	C404-281113-8	13mm	C404-320913-8	13mm	C404-321113-8
14mm	C404-280914-8	14mm	C404-281114-8	14mm	C404-320914-8	14mm	C404-321114-8
15mm	C404-280915-8	15mm	C404-281115-8	15mm	C404-320915-8	15mm	C404-321115-8
16mm	C404-280916-8	16mm	C404-281116-8	16mm	C404-320916-8	16mm	C404-321116-8



For more information or to place an order call: 877.9CURITEVA, email: customersupport@curiteva.com or visit www.curiteva.com

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2020 by Curiteva, Inc.
QMSF-1150 Rev 1.0