


SAVANT™ PEEK

Transforaminal Lumbar Interbody Fusion System
Featuring Titanium Plasma Coating

LUMBAR INTERBODY


Features & Benefits



CURITEVA®

SAVANT™ PEEK TLIF

Featuring Titanium Plasma Coating



The Savant Titanium Plasma Coated Transforaminal Lumbar Interbody Fusion System is designed to offer immediate stability and long-term fixation. The interbody is manufactured out of PEEK with a titanium plasma spray coating. The titanium porous coating creates topographical properties and provides fluoroscopic endplate visualization. The interbody is offered in a variety of footprints, heights and angles.

System Features:

- Surface topography produces superior implant contact and construct stability
- Modulus of elasticity between cortical and cancellous bone for optimal load sharing
- Designed to ensure a high degree of compressive strength and dimensional stability
- Optimal Radiolucency
- Ergonomically shaped anterior edges
- Oversized, central lateral graft windows
- Comprehensive discectomy instruments
- Full range of Trials for proper sizing of implant height
- Tantalum radiopaque markers to optimize visibility and placement
- Multi-directional teeth resist migration
- Self-distracting leading edge to facilitate implant insertion

Titanium Plasma Coated PEEK TLIF Interbodies

28mm Length, 9mm Width, 0°		28mm Length, 11mm Width, 0°		32mm Length, 9mm Width, 0°		32mm Length, 11mm Width, 0°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
6mm	C424-280906-OCT	6mm	C424-281106-OCT	6mm	C424-320906-OCT	6mm	C424-321106-OCT
7mm	C424-280907-OCT	7mm	C424-281107-OCT	7mm	C424-320907-OCT	7mm	C424-321107-OCT
8mm	C424-280908-OCT	8mm	C424-281108-OCT	8mm	C424-320908-OCT	8mm	C424-321108-OCT
9mm	C424-280909-OCT	9mm	C424-281109-OCT	9mm	C424-320909-OCT	9mm	C424-321109-OCT
10mm	C424-280910-OCT	10mm	C424-281110-OCT	10mm	C424-320910-OCT	10mm	C424-321110-OCT
11mm	C424-280911-OCT	11mm	C424-281111-OCT	11mm	C424-320911-OCT	11mm	C424-321111-OCT
12mm	C424-280912-OCT	12mm	C424-281112-OCT	12mm	C424-320912-OCT	12mm	C424-321112-OCT
13mm	C424-280913-OCT	13mm	C424-281113-OCT	13mm	C424-320913-OCT	13mm	C424-321113-OCT
14mm	C424-280914-OCT	14mm	C424-281114-OCT	14mm	C424-320914-OCT	14mm	C424-321114-OCT
15mm	C424-280915-OCT	15mm	C424-281115-OCT	15mm	C424-320915-OCT	15mm	C424-321115-OCT
16mm	C424-280916-OCT	16mm	C424-281116-OCT	16mm	C424-320916-OCT	16mm	C424-321116-OCT

28mm Length, 9mm Width, 8°		28mm Length, 11mm Width, 8°		32mm Length, 9mm Width, 8°		32mm Length, 11mm Width, 8°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
7mm	C424-280907-8CT	7mm	C424-281107-8CT	7mm	C424-320907-8CT	7mm	C424-321107-8CT
8mm	C424-280908-8CT	8mm	C424-281108-8CT	8mm	C424-320908-8CT	8mm	C424-321108-8CT
9mm	C424-280909-8CT	9mm	C424-281109-8CT	9mm	C424-320909-8CT	9mm	C424-321109-8CT
10mm	C424-280910-8CT	10mm	C424-281110-8CT	10mm	C424-320910-8CT	10mm	C424-321110-8CT
11mm	C424-280911-8CT	11mm	C424-281111-8CT	11mm	C424-320911-8CT	11mm	C424-321111-8CT
12mm	C424-280912-8CT	12mm	C424-281112-8CT	12mm	C424-320912-8CT	12mm	C424-321112-8CT
13mm	C424-280913-8CT	13mm	C424-281113-8CT	13mm	C424-320913-8CT	13mm	C424-321113-8CT
14mm	C424-280914-8CT	14mm	C424-281114-8CT	14mm	C424-320914-8CT	14mm	C424-321114-8CT
15mm	C424-280915-8CT	15mm	C424-281115-8CT	15mm	C424-320915-8CT	15mm	C424-321115-8CT
16mm	C424-280916-8CT	16mm	C424-281116-8CT	16mm	C424-320916-8CT	16mm	C424-321116-8CT



For more information or to place an order call: 877.9CURITEVA, email: customersupport@curiteva.com or visit www.curiteva.com

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2020 by Curiteva, Inc.
QMSF-1151 Rev 1.0