


SAVANT™ Titanium

Transforaminal Lumbar Interbody Fusion System
Featuring TRU-LOK® Technology

LUMBAR INTERBODY

Features & Benefits




CURITEVA®

SAVANT™ Titanium TLIF SYSTEM

Featuring TRU-LOK® Technology

Savant Titanium Transforaminal Lumbar Interbody System features TRU-LOK®. The TRU-LOK® surface technology results in a proprietary micro-textured surface intended to provide an improved osteogenic environment. The micro-scale roughness of the TRU-LOK® surface technology is intended to reduce the potential for implant migration and to encourage enhanced osseointegration. The interbody is offered in multiple footprints, heights and lordosis offering an ideal fit for every patient. The Savant Interbody Fusion System features instrumentation that allows for better control and more precise implant placement.



System Features:

- Open Architecture
- Convex design offers a more accurate fit to the patient's anatomy
- Micro-roughness surface texturing produces superior implant contact and construct stability
- Oversized central and lateral graft windows for biologic seeding
- Contoured, distracting nose to facilitate placement
- Multi-directional teeth resist migration
- Comprehensive discectomy instruments
- Full range of Trials for proper sizing of implant height

Titanium TLIF Interbodies with TRU-LOK® Technology

28mm Length, 9mm Width, 0°		28mm Length, 11mm Width, 0°		32mm Length, 9mm Width, 0°		32mm Length, 11mm Width, 0°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
6mm	C414-280906-0	6mm	C414-281106-0	6mm	C414-320906-0	6mm	C414-321106-0
7mm	C414-280907-0	7mm	C414-281107-0	7mm	C414-320907-0	7mm	C414-321107-0
8mm	C414-280908-0	8mm	C414-281108-0	8mm	C414-320908-0	8mm	C414-321108-0
9mm	C414-280909-0	9mm	C414-281109-0	9mm	C414-320909-0	9mm	C414-321109-0
10mm	C414-280910-0	10mm	C414-281110-0	10mm	C414-320910-0	10mm	C414-321110-0
11mm	C414-280911-0	11mm	C414-281111-0	11mm	C414-320911-0	11mm	C414-321111-0
12mm	C414-280912-0	12mm	C414-281112-0	12mm	C414-320912-0	12mm	C414-321112-0
13mm	C414-280913-0	13mm	C414-281113-0	13mm	C414-320913-0	13mm	C414-321113-0
14mm	C414-280914-0	14mm	C414-281114-0	14mm	C414-320914-0	14mm	C414-321114-0
15mm	C414-280915-0	15mm	C414-281115-0	15mm	C414-320915-0	15mm	C414-321115-0
16mm	C414-280916-0	16mm	C414-281116-0	16mm	C414-320916-0	16mm	C414-321116-0

28mm Length, 9mm Width, 8°		28mm Length, 11mm Width, 8°		32mm Length, 9mm Width, 8°		32mm Length, 11mm Width, 8°	
Height	Catalog Number	Height	Catalog Number	Height	Catalog Number	Height	Catalog Number
7mm	C414-280907-8	7mm	C414-281107-8	7mm	C414-320907-8	7mm	C414-321107-8
8mm	C414-280908-8	8mm	C414-281108-8	8mm	C414-320908-8	8mm	C414-321108-8
9mm	C414-280909-8	9mm	C414-281109-8	9mm	C414-320909-8	9mm	C414-321109-8
10mm	C414-280910-8	10mm	C414-281110-8	10mm	C414-320910-8	10mm	C414-321110-8
11mm	C414-280911-8	11mm	C414-281111-8	11mm	C414-320911-8	11mm	C414-321111-8
12mm	C414-280912-8	12mm	C414-281112-8	12mm	C414-320912-8	12mm	C414-321112-8
13mm	C414-280913-8	13mm	C414-281113-8	13mm	C414-320913-8	13mm	C414-321113-8
14mm	C414-280914-8	14mm	C414-281114-8	14mm	C414-320914-8	14mm	C414-321114-8
15mm	C414-280915-8	15mm	C414-281115-8	15mm	C414-320915-8	15mm	C414-321115-8
16mm	C414-280916-8	16mm	C414-281116-8	16mm	C414-320916-8	16mm	C414-321116-8



For more information or to place an order call: 877.9CURITEVA, email: customersupport@curiteva.com or visit www.curiteva.com

Products patented and/or patent pending.
All products are not currently available in all markets.

Copyright©2020 by Curiteva, Inc.
QMSF-1152 Rev 1.0